TIBBİ CİHAZ / MEDICAL DEVICE

SAÇ EKİM PENSETLERİ

HAIR TRANSPLANT FORCEPS

KULLANMA KILAVUZU / INSTRUCTIONS FOR USE

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01 Product description and intended use

Hair transplant forceps are precision manually operated surgical instruments intended to assist with the collection, arrangement, handling and placement of grafts (follicular units) during hair transplantation procedures.

The products are manufactured from 301 and 304 stainless steel. They are supplied non-sterile and must be cleaned and sterilized before use.

FORCEPS TYPE	INTENDED USE
Gripping forceps	Collection of grafts/hair follicles during hair transplantation.
Implantation forceps	Holding and placing harvested grafts during the implantation stage.
Arrangement forceps	Arranging grafts within a Petri dish.
Curved forceps	Collection of grafts/hair follicles during hair transplantation.

02 Intended users and conditions of use

The devices must be used by authorized healthcare professionals trained in hair transplantation and safe use of medical instruments. Determination of the patient's suitability for hair transplantation is the responsibility of the treating healthcare professional.

IMPORTANT SAFETY INFORMATION

Do not use forceps that are damaged, corroded, bent, misaligned, functionally defective, or have not been adequately cleaned and sterilized.

03 Inspection before use

- Select the appropriate forceps type and size; check the product label and traceability information.
- Inspect the tips, surfaces and gripping function for deformation, cracks, wear, stains or corrosion.
- Remove newly purchased forceps from the packaging and clean and sterilize them before first use.
- Prepare the sterilized instrument using aseptic technique. Never use a device with incomplete sterilization.

04 Precautions during use

- Grip each graft using the minimum force needed for the task and in accordance with the tip design.
- Avoid excessive pressure and force that could crush or traumatize grafts.
- Do not drop the forceps, strike them against hard surfaces or expose their delicate tips to mechanical stress.
- Do not use the device for cutting, piercing, levering or other purposes outside the intended use.
- Remove the instrument from service if damage or impaired function is observed during use.

EN / INSTRUCTIONS FOR USE / REPROCESSING & TRACEABILITY

05 Cleaning and maintenance after use

- Transport used instruments with due regard to contamination risk and follow the institution's infection prevention procedures.
- After the procedure, remove soil and biological residues with suitable cleaning agents.
- Avoid chloride-containing cleaning agents, corrosive chemicals and abrasive metal brushes or instruments.
- Use demineralized (purified) water for cleaning and rinsing; rinse thoroughly and dry the instrument.
- After cleaning, inspect the surfaces, tips and function again and segregate defective instruments.

06 Sterilization

The products are supplied non-sterile. After completion of cleaning, they must be steam sterilized in an autoclave before first use and before each subsequent reuse.

STERILIZATION METHOD	TEMPERATURE	TIME
Steam autoclave	121 °C	30 minutes

Use a suitable sterilization basket, silicone support or container to protect the delicate tips. Arrange the instruments to prevent impacts against hard surfaces during packaging, loading and handling.

07 Reuse and service life

The products are reusable. Provided that they are regularly maintained and used according to their intended purpose, the recommended service life is 2 (two) years. Repeated-use studies cover 105 uses of the product.

Before each use, check integrity, corrosion, tip alignment and gripping function. Record the number of uses for each instrument (see the instrument use record). Remove the instrument from service if damage or functional deterioration is observed, regardless of its period of service.

08 Storage, transport and disposal

Store the forceps in a clean, dry, adequately ventilated area protected from dust, pests, direct sunlight and corrosive vapors. Transport used devices with due regard to contamination and injury risk. Dispose of unserviceable devices according to the healthcare facility's medical waste and sharps procedures.

09 Label information and incident reporting

Observe the symbols that may appear on the product label, including MD (medical device), REF (reference number), LOT (batch number), UDI (unique device identifier), manufacturer, date of manufacture, NON STERILE and consult instructions for use.

Report any serious incident related to the device to the manufacturer and the relevant competent authority. Record the REF, LOT and, where available, UDI information for incident investigation.

10 Manufacturer

ER TIBBİ ÜRÜNLER SAĞLIK HİZ. PAZ. VE DIŞ. TİC. LTD. ŞTİ.
Halide Edip Adıvar Mah. Akar Cad. No:38 D:8 Kat:3, Bomonti - Şişli / İstanbul, Türkiye
Phone: +90 212 210 32 20 | E-mail: info@ertip.com | SRN: TR-MF-000027136

TR + EN / KULLANIM TAKİP FORMU / INSTRUMENT USE RECORD

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Instrument use and reprocessing count

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